



FURIX

FUP-BW-2A

**Portable Veterinary Machine
Ultrasound System**



Product Features

01 15-inch full-view medical LCD display

- Dual probe interface allows simultaneous use of two probes
- Built-in lithium battery, power off for 5 hours
- Adopt ARM chip architecture to provide stable and simple operating system
- The host automatically recognizes and uses multiple array probes
- Support harmonic imaging
- B, B+B, B+M, M, 4B display modes
- Special software for small animals - supports routine measurement of baby date, heart rate and other tissues and organs of small animals
- Graphic management system - check images can be directly stored, exported, browsed, automatically import selected diagnostic images and check measurement data
- Languages: Chinese, English, French, Spanish, German, Portuguese, Russian, Turkish, Arabic, customized languages according to different users
- Combined with a large number of clinical experience, the machine
- has different preset adjustment conditions for different parts or organs examined by different departments and has a one-click optimization function



3.5MHz Convex Probe





5MHz Micro-convex probe



7.5MHz Linear Probe



Trans-rectal probe



Specification

Monitor:	15 inch LCD medical monitor
Platform:	80
RAM:	128GB
Probe connector:	2
RAM:	128GB
Port:	S-Video, USB, Video, VGA, LAN, Foot Switch
Language:	Chinese, English, Spanish, French, German, Czech, Russian etc
Application:	Pet abdomen, heart, obstetrics, small organs, blood vessels, muscles and bones etc



Shenzhen Falcon and Furix Lifesciences Co. Ltd
Furix Trade Venture Limited (HK)

Address: 6/F MEGA CUBE, WANG KWONG ROAD, KOWLOON, HONG KONG,
 Mob:+852 9579 0472 E-mail: vet.furix@gmail.com, sales@myfurix.com
 Tel:+852 9579 0472

FURIX