



FURIX

FC-05MC

FURIX MEDICAL CENTRIFUGE

Intelligent control, low noise and stability Small size, large capacity, small footprint



Medical Centrifuge

- Maintenance-free brushless motor, LCD display speed and time
- Adjustable speed, from 300 to 5000rpm
- Timer from 155 to 99min or continuous
- 4×50ml/15ml or 4×1.5/2.0ml Suitable for 4×10ml /7ml/5ml/3ml blood collection tube, 4×50ml/15ml centrifugal tube or 4×1/5/2.0ml
- The lid opens automatically at the end of the run. To prevent sample warming and shorten processing time
- Two program buttons P1 / P2,easily access to regular programs
- Easily switch between RPM and RCF with the press of a button



Pet Healthcare

Technical Parameter

MODEL	FC-05MC
SPEED RANGE	300-5000RPM
MAX RCF	2500XG
TIME	15S-99MIN OR CONTINUOUS
POWER SUPPLY	100-240V 50/60HZ / DC24V
CAPACITY	4X50ML P2 15ML/10ML/7ML/5ML/3ML/1.5ML 6X20ML OPTIONAL / 8X10ML OPTIONAL
DIMENSION (MM)	295X260X185
PACKAGE (MM)	350X310X250
G.W./N.W.	6.5 / 5.5KG

Furix

Application



Widely used in scientific research, education and production departments such as biology, chemistry, medicine, etc., suitable for hospitals, blood banks, clinical trials and biochemical laboratories and other institutions



Furix Trade Venture Limited

Address: 6/F MEGA CUBE, WANG KWONG ROAD, KOWLOON, HONG KONG,
Mob: +852 9579 0472 E-mail: vet.furix@gmail.com, sales@myfurix.com
Tel: +852 9579 0472 www.myfurix.com

FURIX