



FURIX

FWF

Fixed Wooden Frame
For Pet Treatment



Product Features

- 01 Pet specific protective frame, including three major parts: main frame, bilateral chest rest, and protective strap
- 02 There are knobs on both sides of the chest support, with adjustable width to adapt to the patient's chest width
- 03 The main bracket is made of solid wood, wrapped in a leather cover and lining, providing comfortable sleeping experience
- 04 The fixed strap can adjust the tightness and length, suitable for patients of different sizes
- 05 The animals conduct safe acupuncture and moxibustion operations in a conscious state, fully exposing commonly used acupuncture and moxibustion points
- 06 Suitable for other operations such as blood collection, injection, infusion, wound cleaning, local anesthesia surgery, etc

Parameter



Standard model	Length*width*height	reference variety	reference weight
S	50*28*42	Baby pets, small dogs, rabbits, cats	1-10 KG
M	60*30*48	Small and medium-sized dogs (such as Chai Dog and Silver Fox)	5-20 KG
L	75*33*65	Medium dog, large dog (such as golden coat, Labrador)	15-35 KG



Furix Trade Venture Limited

Address: 6/F MEGA CUBE, WANG KWONG ROAD, KOWLOON, HONG KONG,
Mob:+852 9579 0472 E-mail: vet.furix@gmail.com, sales@myfurix.com
Tel:+852 9579 0472 www.myfurix.com

FURIX