



FURIX

FOC-10W

Dual flow

Furix Oxygen Concentrator

DAILY OXYGEN PROTECTION AND HEALTH CARE



FOC-10W

- **Advanced Humidifier:** Features a USA-imported humidifier with 6PSI pressure release for enhanced performance.
- **Big LCD Display:** Displays key information such as switch times, operating pressure, present working time, accumulated working time, and pre-setting time from 10 minutes to 40 hours.
- **Continuous Operation:** Designed for 24/7 operation throughout the year, ensuring consistent oxygen supply.
- **Innovative Cooling Technology:** Enhances reliability and extends the operational life of the unit.
- **Wide Flow Rate Range:** Supports a flow rate from 1 to 10 LPM, catering to varying oxygen needs.
- **High Purity Oxygen:** Delivers oxygen with 93% purity ($\pm 3\%$), suitable for medical-grade requirements.
- **Low Noise Operation:** Operates at a sound level of ≤ 50 dB for a quiet and comfortable experience.
- **Built-in Alarm System:** Includes power failure, high-pressure, and low-pressure alarms for safety and reliability.



Pet Healthcare Oxygen Concentrator

Specification


Flow Rate	1-10LPM
Purity	93% ($\pm 3\%$)
Outlet Pressure	0.04-0.07MPa
Sound Level	≤ 50 db
Power Consumption	≤ 880 w
LCD Display	Switch times, operating pressure, present working time, accumulated working time, presetting time from 10mins to 40hours
Product Weight	27 Kgs
Product Size	365×375×600mm
Gross Weight	29 Kgs
Carton Size	435×435×640mm



USA imported humidifier
6PSI pressure release



Innovative "cooling"
technology to improve
reliability and life of
the unit



Big LCD display
24H*365D
work continuously

Imported molecular sieve gives you a full sense of security with natural enjoyment!

Furix

— Furix Oxygen Concentrator FOC-10W



Single flow

Furix Trade Venture Limited

Address: 6/F MEGA CUBE, WANG KWONG ROAD, KOWLOON, HONG KONG,
Mob: +852 9579 0472 E-mail: vet.furix@gmail.com, sales@myfurix.com
Tel: +852 9579 0472 www.myfurix.com

FURIX